



## Detecto 6856 High Capacity Digital Handrail Scale, 1,000 lb x 0.2 lb

Weighs in lbs & kgs, with BMI feature and RS232 interface, includes new Detecto 758C indicator, large 24 x 24 inch low profile platform, auto zero, display hold, motion detection. Battery or optional AC adapter



Manufacture: Detecto  
SKU: 6856  
UPC: 809161137408

**Please Call  
(800)832-0055**

### Features

Heavy-duty, yet lightweight, the 6856 bariatric scales balanced design makes it easy to move. The 6856 is engineered specifically for weighing obese and unsteady patients and features a 24 (610mm) square, low-profile platform with non-skid surface and sturdy handrails with padded grips.

The digital indicator is the industrys most advanced, with easy-to-use quick weigh keys, plus Body Mass Index calculation. The auto-zero function ensures accurate weight each time without having to press the zero key; the display hold feature retains weight so you can attend to the patient first, then record the weight. Powered by your choice of either batteries or AC adapter, the 6856 also has an auto-off function to conserve energy.

#### Specifications:

- Capacity: 800 lbs x 0.2 lbs (360 kg x 0.1 kg)
- Includes NEW Detecto 750 indicator
- Highly accurate digital weight display
- Body Mass Index feature
- RS-232 interface
- Extra strength wrap-around tubular stainless steel handrails
- LCD weight display which can be read from either side of the platform
- Large heavy-duty rubber wheels and the scale tips back for easy transporting
- Skid-proof and soft, but durable black rubber mat
- Platform Size: 24 x 24 inches (610 x 610 mm)
- Power: AC adapter or 6 AA size alkaline batteries
- Display: 0.7 / 18 mm high-contrast, 7-segment LCD with Body Mass Index, On/Off button, Lb/Kg, Hold/Release, Zero, Automatic Zero
- Display Functions: Auto power off, Hold/release of display, Motion detection

Power x 3 with the new 750 Weight Display

- The new Detecto 750 weighing indicator is the perfect weight display for all applications. The recharging circuitry is “Built-In” to the 750 display. Utilizing the AC-adapter and installing Ni-Cad “AA” batteries, you can have a truly free standing scale in any office environment.
- Operate up to 250 hours on new alkaline batteries.
- Operate up to 100 hours on NiCad batteries before recharging. (Recharging time is approximately 8 hours and the scale can be operated during this time.)

## Accessories

### Optional Accessories

---



### Detecto AC Adapter

AC Adapter for Detecto Medical Scales

---

[Contact US](#)

[718-336-5900](#)

[Returns](#)

[Store Policies](#)

[About Us](#)

- Phone: [\(800\)832-0055](#) or [\(718\)336-5900](#)
- Chat with us via [WhatsApp Chat](#)
- Address: 4802 Glenwood Rd. Brooklyn, NY 11234
- © 2000 - 2020 scalesgalore.com (a division of Itin Scale Co., Inc.) All rights reserved.