

Detecto's P50 printer can be mounted or used with any ProDoc medical scales that have RS232 serial connectivity. AC power and direct thermal dot line printing make the P50 ideal for printing patient records.



**Manufacture:** Detecto

**SKU:** P50

**Dimensions**

**X:** 4.33 **Y:** 4.33 **Z:** 4.21

**Please Call**  
**(800)832-0055**

## Features

### P50 Direct Thermal Printer With Serial Interface

Detectos P50 printer can be mounted or used with any ProDoc medical scales that have RS232 serial connectivity. AC power and direct thermal dot line printing make the P50 ideal for printing patient records.

Print Time, Date, Height, Weight, Age, Gender, Activity Level, BMI, Body Fat, Body Water, Bone Mass, and Muscle Mass on DETECTOs ProDoc Body Fat Analysis Scales.

## Specifications

P50 Specifications	
<b>Printing Type:</b>	Direct Thermal Dot Line
<b>Number of Dots/Line:</b>	384 dots
<b>Resolution:</b>	203 dpi (8 dots/mm W x 16 dots/mm H)

i81u812

Contact Us 718-336-5000 Returns  
 Store Policies About Us

<b>Max. Printing Speed:</b>	1.58 in/second (40 mm/second)
<b>Max. Paper Roll Diameter:</b>	1.97 in (50 mm)
<b>Printing Width:</b>	1.89 in (48 mm)
<b>Operating Temperature:</b>	68 to 122° F (20 to 50° C)
<b>Power Supply:</b>	AC Adapter 100 - 240VAC input, 12V [2A] output
<b>Connectivity:</b>	RS-232 communication with scale (cable included)
<b>Dimensions:</b>	4.33 in W x 4.33 in D x 4.21 in H (110 mm x 110 mm x 107 mm)

- 
- Phone: [\(800\)832-0055](tel:8008320055) or [\(718\)336-5900](tel:7183365900)
  - Chat with us via [WhatsApp Chat](#)
  - Address: 4802 Glenwood Rd. Brooklyn, NY 11234
  - © 2000 - 2024 scalesgalore.com (a division of Itin Scale Co., Inc.) All rights reserved.

P50 Features:
• Direct thermal printing
• Compact design
• RS-232 serial connectivity - includes compact cable ready to plug into Detectos ProDoc scales
• Two-button operation
• Thermal dot line printing
• Easy to change paper roll
• Wall mounting or may be attached to ProDoc eye-level scales
• Error and power indicators
Serial Connectivity
Use it with the ProDoc Series Scales
Use it with the ProDoc Body Fat Analysis Series Scales