

Suitable for normal samples. Use with OHAUS pH meters with BNC input connector and Cinch temperature connector.



Manufacture: Ohaus

SKU: ST310

Weight: 1.00 lb

Free Ground Shipping
within the 48
continental US States

Please Call
(800)832-0055

Features

3 in 1 ST310 Plastic Refillable pH Electrode

This delicate measuring instrument requires proper care and maintenance to produce accurate and reliable results and to prolong useful life. Always keep the electrode moist when not in use for a period of time by using an electrode storage solution (3M KCl). DO NOT store the electrode in distilled or deionized water, as this will cause ions to leach out of the glass bulb and reference electrolyte, causing slow and sluggish response. Remove the electrode gently from the storage bottle and rinse it with distilled water before use. For long-term storage, always keep the electrode in the bottle filled with sufficient storage solution to cover the bulb. Replenish the bottle as needed.

Suitable for normal samples. Use with OHAUS pH meters with BNC input connector and Cinch temperature connector.

	ST310	ST210	ST320
pH Range	0 to 14 pH	0 to 14 pH	0 to 14 pH
Temp Range	0 to 80°C	0 to 80°C	0 to 80°C
Shaft Material	Plastic	Plastic	Plastic
Internal Reference Type	Ag/AgCl	Ag/AgCl	Ag/AgCl
Gel/Refillable	Refillable	Refillable	Non-Refillable Gel

Reference Junction Type	Ceramic pin	Ceramic pin	Fiber pin
Refilling Reference Electrolyte	3M KCl solution	3M KCl solution	3M KCl gel
Shaft Dimensions	4-3/4" x 3/4"	4-3/4" x 3/4"	4-3/4" x 3/4"
Cable Length	3-1/4"	3-1/4"	3-1/4"
Temp Sensor	Yes	No	Yes
Connector	BNC Cinch	BNC	BNC

i81u812

Contact US

718-336-5900

Returns

Store Policies

About Us

- Phone: [\(800\)832-0055](tel:8008320055) or [\(718\)336-5900](tel:7183365900)
 - Chat with us via [WhatsApp Chat](#)
- Address: 4802 Glenwood Rd. Brooklyn, NY 11234
- © 2000 - 2024 scalesgalore.com (a division of Itin Scale Co., Inc.) All rights reserved.